

NEW GENERATION MESHED MAT



injection



surgery



orthopedics



hospitalization



capture



protection



cleaning



livestock



grooming


génia
the veterinarian's preference

NEW

COSYPAD®

New generation meshed mat realized in a material stemming from the medical environment used as anti-bedsores, Cosypad® offers comfort, hygiene and ease of cleaning. It constitutes an alternative of choice in duck-boards, underpads, and others litters for the well-being of animals and better efficiency on a daily basis. Thanks to its aerated structure, animals always stay clean and dry. Cosypad® satisfies all the conditions allying the comfort for our animals and the gain of time for professionals.

- Reduces the pressure exercised on muscles and articulations for an optimal comfort and the prevention of bedsores.
- Hygienic, liquids are evacuated under the mat leaving the animal dry.
- Do not retain hairs.
- Good aeration allowing a better healing of wounds and a better thermoregulation.
- Hypoallergenic material, odorless, antibacteria.
- Very resistant, rot-proof and insulating.
- Reinforced stitching. Super resistant nylon on the edges for better durability.
- One non-slip surface.
- Cleaning instructions:
 - in a washing machine at 40°C with a bleach tablet
 - with a disinfection spray by soaking in an appropriate solution
 - spinning/dryer not recommended
- One suspension eye to simplify storage.
- Uses:
 - in hospitalization kennels
 - on consultation/surgery tables
 - during post-operative wake up stage
 - for descaling (fast evacuation of liquids)
 - as anti-fatigue mat (for the surgeon)



Code	Size (cm)	Weight (g)
327200	33 x 50	200
327201	39 x 71	382
327203	55 x 71	505
327205	60 x 115	884



COMFORT, HYGIENE | GOOD AERATION | STURDY & DURABLE MATERIAL | & MACHINE WASHABLE

TESTIMONIES

Numerous veterinary practices already use the concept and show innovative qualities of the pad.

Alfort Veterinary School (France)

- Excellent thermal insulation to use instead of duck-boards and underpads in cats kennels.
- A better grip for animals, thus more safety especially after orthopedic surgeries.
- More hygiene because animals are not soiled by urines anymore.
- Provide an optimal comfort in dogs kennels.
- Prevent the formation of bedsores on animals during long hospitalization.
- Easy washing at 40°C with a disinfectant

Armand Valla – Rauzan (France)

We chose to use Cosypad® in surgery room because their characteristics are optimal for the following needs : isolation of animals on the table of surgery, management of the hypothermia, perfect fluids flow under the animal which remains dry and not soiled, comfort and stability of the animals placed in dorsal decubitus. Pads are daily washed in machine with bleach tablets, dry extremely fast and present no sign of wear.

Dr. Jean-paul Zendali – Fribourg (Switzerland)

I use Cosypad® in the recovery rooms, after surgeries. Cats or dogs always urinate in phase of awakening – because of pre-surgery perfusion - and it is always unpleasant to return an animal still wet to its owner. With Cosypad®, this problem is solved ! The cleaning is easy : under the shower or in washing machine. It is comfortable and soft for animals in the awakening, which never damaged it so far! Really an advantage for the daily use.



HOSPITALIZATION